

# Arraystar G4 CUT&Tag Library Prep Kit

Cat#: AS-TN-006

## Instruction Manual version 1.0

For research use only. Not for diagnostic procedures.

### Background

G-quadruplexes (G4s) are four-stranded intramolecular structures that arise from the self-stacking of two or more guanine quartets (G-quartets). G4s are distributed in specific genomic regions and implicated in several essential cellular processes. This kit uses G4 CUT&Tag (Cleavage Under Targets and Tagmentation) to profile native DNA G4s at high resolution, specificity, and sensitivity.

### Components

| Component                    | Volume    | Storage | Bag   |
|------------------------------|-----------|---------|-------|
| Wash Buffer                  | 1.5 mL    | 4°C     | Bag 1 |
| Antibody Buffer              | 1.5 mL    |         |       |
| Dig-150 Buffer (10x)         | 2.6 mL    |         |       |
| Dig-300 Buffer (10x)         | 1.5 mL    |         |       |
| Tagmentation Buffer          | 20 uL     |         |       |
| Stop Buffer                  | 80 uL     |         |       |
| Neutralize Buffer            | 20 uL     |         |       |
| 60 mM Spermidine             | 400 uL    |         |       |
| 1% NP40                      | 500 uL    |         |       |
| DNA Beads                    | 2 mL      |         |       |
| Antibody 1                   | 75 uL     | -80°C   | Bag 2 |
| Assembled pA-Tn5 Transposome | 25 uL     | -20°C   |       |
| Antibody 2                   | 14 uL     |         |       |
| Antibody 3                   | 14 uL     |         |       |
| P5/P7 PCR Primers (20 uM)    | 6 x 12 uL |         |       |
| PCR Mix (2x)                 | 280 uL    |         |       |
| 1% Digitonin                 | 2 mL      |         |       |

**Note:** Store Antibody 1 separately at -80°C upon receipt. Avoid freeze-thaw cycles for 1% Digitonin.

### Additional Materials

- Nuclease-free ddH<sub>2</sub>O
- 80% ethanol
- BSA dry powder
- Vortexer
- Magnetic plate rack
- Centrifuge
- Orbital shaker
- Water bath or incubator
- Thermal cycler
- Analytical balance

### Important Notes

- This kit is suitable for cultured cell or cells isolated from tissue samples. The recommended cell viability is > 85%.
- To avoid cell loss during wash processes, leave 3~5 uL liquid when removing the supernatant. A horizontal centrifuge is easier for cell collection.

### G4 CUT&Tag Protocol

The entire protocol takes three days to complete.

#### DAY 1

Prepare fresh buffers according to the tables below. Keep the buffers on ice.

| 10% BSA                          |                   |
|----------------------------------|-------------------|
| Component                        | Amount per sample |
| BSA dry power                    | 50 mg             |
| Nuclease-free ddH <sub>2</sub> O | 500 uL            |

**Note:** Store the freshly prepared BSA at 4°C. It will be used for buffer prep on Day 1 and Day 2.

| Complete Wash Buffer |                   |
|----------------------|-------------------|
| Component            | Amount per sample |
| Wash Buffer          | 100 uL            |
| 60 mM Spermidine     | 0.8 uL            |

| Complete Antibody Buffer |                   |
|--------------------------|-------------------|
| Component                | Amount per sample |
| Antibody Buffer          | 83 uL             |
| 1% Digitonin             | 5 uL              |
| 1 % NP40                 | 1 uL              |
| 60 mM Spermidine         | 0.8 uL            |
| 10% BSA                  | 10 uL             |

## 1. Incubate Antibody 1

1.1. Count the cells and check the cell viability (> 85%) by Trypan Blue staining. Aliquot 200,000 cells to a 1.5 mL tube per sample.

1.2. Centrifuge the cells at 600 x g for 3 minutes at 4°C. Remove and discard the supernatant.

1.3. Resuspend the cells with 100 uL Complete Wash Buffer. Centrifuge cells at 600 x g for 3 minutes at 4°C. Remove and discard the supernatant.

1.4. Resuspend the cells with 94 uL of Complete Antibody Buffer and 6 uL of Antibody 1. Incubate the tube on an orbital shaker at 20 RPM overnight at 4°C.

## DAY 2

Prepare fresh buffers according to the tables below. Keep the buffers on ice.

| Complete Dig-150 Buffer          |                   |
|----------------------------------|-------------------|
| Component                        | Amount per sample |
| Dig-150 Buffer (10x)             | 210 uL            |
| 1% Digitonin                     | 105 uL            |
| 60 mM Spermidine                 | 17.5 uL           |
| 10% BSA (from Day 1)             | 210 uL            |
| 1 % NP40                         | 21 uL             |
| Nuclease-free ddH <sub>2</sub> O | 1537 uL           |

| Complete Dig-300 Buffer          |                   |
|----------------------------------|-------------------|
| Component                        | Amount per sample |
| Dig-300 Buffer (10x)             | 110 uL            |
| 1% Digitonin                     | 55 uL             |
| 60 mM Spermidine                 | 9.2 uL            |
| 1 % NP40                         | 11 uL             |
| 10% BSA (from Day 1)             | 110 uL            |
| Nuclease-free ddH <sub>2</sub> O | 805 uL            |

## 2. Incubate Antibody 2

2.1. Centrifuge cells at 600 x g for 5 minutes at 4°C. Remove and discard the supernatant. Avoid disturbing the cells.

2.2. Resuspend the cells with 200 uL Complete Dig-150 Buffer by gentle pipetting. Incubate for 10 minutes at room temperature, gently mixing the cells once with a pipette at the 5-minute mark.

2.3. Repeat the wash steps 2.1 - 2.2 twice (three washes in total).

2.4. Resuspend the cells with 99 uL Complete Dig-150 Buffer and 1 uL Antibody 2 by gentle pipetting. Transfer to a new 1.5 mL tube. Incubate the tube on an orbital shaker at 20 RPM for 1 hour at room temperature.

## 3. Incubate Antibody 3

3.1. Repeat wash steps 2.1 - 2.2 three times in total.

3.2. Resuspend the cells with 99 uL Complete Dig-150 Buffer and 1 uL Antibody 3 by gentle pipetting. Transfer to a new 1.5 mL tube. Incubate the tube on an orbital shaker at 20 RPM for 1 hour at room temperature.

## 4. Incubate pA-Tn5 transposome

4.1. Repeat wash steps 2.1 - 2.2 three times in total.

4.2. Resuspend the cells with 98 uL Complete Dig-300 Buffer and 2 uL Assembled pA-Tn5 Transposome by gentle pipetting. Transfer to a new 1.5 mL tube. Incubate the tube on an orbital shaker at 20 RPM for 1 hour at room temperature.

## 5. Tagmentation of target DNA

5.1. Centrifuge the cells at 600 x g for 5 minutes at 4°C. Remove and discard the supernatant. Avoid disturbing the cells.

5.2. Resuspend the cells with 300 uL Complete Dig-300 Buffer by gentle pipetting. Incubate for 10 minutes at room temperature, gently mixing the cells once with a pipette at the 5-minute mark.

5.3. Repeat wash steps 5.1 - 5.2 twice (three washes in total).

5.4. Resuspend the cells with 49.5 uL Complete Dig-300 Buffer and 0.5 uL Tagmentation Buffer by gentle pipetting. Transfer to a new 1.5 mL tube. Incubate at 37°C for 1 hour, with gentle mixing of the cells once using a pipette at the 30-minute time point.

## 6. Release of DNA

6.1. Mix Stop Buffer just before use. Add 6 uL Stop Buffer and 1 uL of Neutralize Buffer by gentle pipetting. Incubate at 55°C for 1 hour, gently mixing the cells with a pipette at the 20- and 40-minute time points.

6.2. Add 110 uL DNA Beads and vortex to mix. Incubate at room temperature for 10 minutes.

6.3. Separate the beads by placing the tube on a magnetic rack for 4 minutes. Remove and discard the supernatant.

6.4. Keep the tube on the magnetic rack, add 350 µL of 80% ethanol without disturbing the beads. After 1 minute, remove and discard the supernatant. Repeat this step once for a total of two washes.

6.5. Air dry at 37°C for 3 min to completely remove the residual ethanol. Resuspend the beads with 20 uL Nuclease-free ddH<sub>2</sub>O and transfer the bead suspension to a new PCR tube.

## 7. Library amplification

7.1. Add 1 uL P5 PCR Primer (20 uM), 1 uL P7 PCR Primer (20 uM), and 22 uL PCR Mix (2x) to the PCR tube by gentle pipetting.

7.2. PCR Amplification using the following PCR program:

| Cycle number         | Temperature | Time       |
|----------------------|-------------|------------|
| 1X (initial heating) | 72°C        | 5 minutes  |
| 1X (denaturation)    | 98°C        | 30 seconds |
| 14 ~ 17X             | 98°C        | 10 seconds |
|                      | 60°C        | 30 seconds |
|                      | 72°C        | 60 seconds |
| 1X                   | 72°C        | 5 minutes  |
| 1X                   | 4°C         | Hold       |

## DAY 3

### 8. Enrichment of library fragments

8.1. Resuspend and transfer the beads to a new 1.5 mL tube. Add an extra 27 uL fresh DNA Beads to PCR product. Mix thoroughly and incubate at room temperature for 10

minutes.

8.2. Place the tube on a magnetic rack for 4 minutes. Transfer the supernatant to a new 1.5 mL tube.

8.3. Add 17 uL fresh DNA Beads to the isolated supernatant. Mix and incubate at room temperature for 10 minutes.

8.4. Place the tube on a magnetic rack for 4 minutes. Remove and discard the supernatant.

8.5. Keep the tube on the magnetic rack, add 350 µL of 80% ethanol without disturbing the beads. After 1 minute, remove and discard the supernatant. Repeat this step once for a total of two washes.

8.6. Allow the beads to air dry completely. Add 10-20 uL Nuclease-free ddH<sub>2</sub>O to resuspend the beads and elute the DNA.

8.7. Use the magnetic rack to separate the beads. Transfer the supernatant (library) to a fresh tube.


## 9. Library quantification

9.1. Take 1 uL library for DNA concentration measurement with a Qubit™ Fluorometer.


- If the total DNA amount is ≥ 20 ng, the library is ready for sequencing.
- If the total DNA amount is < 20 ng, PCR amplify the library with 2-4 more cycles in step 7.

## Kit performance

A G4 CUT&Tag Library for 200,000 HeLa cells was prepared with this kit according to the kit protocol. The library was sequenced on an Illumina Sequencer. The typical distribution of G4 peaks in the genome is displayed in Fig. 1. G4 peaks in CUT&Tag-seq and ChIP-seq correlate well as shown in Fig.2.



**Fig. 1.** The distribution of G4 peaks in the genome, showing characteristic enrichment at promoter.



**Fig. 2.** The IGV view of G4 ChIP (upper) and G4 CUT&Tag (lower) tracks

## Troubleshooting

| Problems                        | Possible reasons                   | Solutions                                    |
|---------------------------------|------------------------------------|--|
| Low signal-to-background ratios | Often caused by unhealthy cells or | 1. Pre-treat cells with DNase or sort viable |

|                           |  |  |
|---------------------------|--|--|
| in low-depth sequencing   | other non-ideal input material.  | cells.<br>2. Deeper sequencing.  |
| Low library amplification | 1. Special sequence or structure difficult for PCR.<br>2. Residual ethanol in Step 6.3 that inhibited PCR. | 1. Increase PCR cycles.<br>2. Add Nuclease-free ddH <sub>2</sub> O to elute the DNA after the ethanol has completely evaporated. |

## Appendixes

### 1. P5 and P7 PCR primer sequences (5'→3')

|             |  |
|-------------|--|
| P5, Index 1 | AATGATACGGCGACCACCGAGATCTACAC TATCCTCT<br>TCGTCGGCAGCGTC |
| P5, Index 2 | AATGATACGGCGACCACCGAGATCTACAC GTAAGGA<br>TCGTCGGCAGCGTC  |
| P5, Index 3 | AATGATACGGCGACCACCGAGATCTACAC ACTGCATA<br>TCGTCGGCAGCGTC |
| P7, Index 1 | CAAGCAGAAGACGGCATACGAGAT CTAGTACG GTCT<br>CGTGGGCTCGG    |
| P7, Index 2 | CAAGCAGAAGACGGCATACGAGAT TCTGCCT GTCTC<br>GTGGGCTCGG     |
| P7, Index 3 | CAAGCAGAAGACGGCATACGAGAT GCTCAGGA GTCT<br>CGTGGGCTCGG    |

**Note:** The sample multiplexing indexes are yellow highlighted and marked in green and red. A P5 can be paired with a P7 PCR primer in any index combination for sample multiplexing in sequencer.

### 2. Library fragment schematic sequence

5'-  
AATGATACGGCGACCACCGAGATCTACAC TATCCTCT  
AGATGTGTATAAGAGACAG NNNNNN  
CTGTCTCTTATACATCTCCGAGCCCACGAGAC TCTCCTATG  
CCGTCTTCTGCTT-3'

TATCCTCT: P5 Index (8 bases)  
TCTCCTATG: P7 Index (8 bases)  
NNNNNN: Target DNA sequence