

Arraystar R-loop CUT&Tag Library Prep Kit

Cat#: AS-TN-005

Instruction Manual version 1.0

For research use only. Not for diagnostic procedures.

Background

R-loops are special RNA:DNA:DNA three-stranded structures that widely exist in the genome. R-loops participate in many biological processes and play important roles in regulating genome stability. Traditional R-loop detection by DRIP-seq is based on immunoprecipitation of ultrasound or enzyme fragmented genomic DNA containing R-loops. Greatly improved upon traditional methods, Arraystar R-loop CUT&Tag Library Prep Kit uses pre-assembled transposomes containing hyperactive Tn5 transposase fusion with protein-A and transposons to target S9.6 antibody bound R-loops and tagment the R-loop sites (CUT&Tag, Cleavage Under Targets and Tagmentation)(Fig. 1). This powerful method vastly simplifies the R-loop sequencing library prep, allowing low sample inputs, easy lab procedures, fast and better results.




Fig.1. Arraystar R-loop CUT&Tag Library Prep Kit uses transposase/protein A fusion to target and tagment R-loops via S9.6 antibody.

Components

Component	Volume	Storage	Bag
Wash Buffer	1.5 mL	4°C	Bag 1
Antibody Buffer	1.5 mL		
Dig-150 Buffer (10x)	2 mL		
Dig-300 Buffer (10x)	1.5 mL		
Tagmentation Buffer	20 uL		
Stop Buffer	80 uL		
Neutralize Buffer	20 uL		
60 mM Spermidine	300 uL		
1% NP40	400 uL		
DNA Beads	2 mL		
Antibody 1	14 uL	-80°C	Bag 2
Assembled pA-Tn5 Transposome	25 uL	-20°C	
Antibody 2	14 uL		
P5/P7 PCR Primers (20 uM)	6 x 12 uL		
PCR Mix (2x)	280 uL		
1% Digitonin	1.8 mL		

Note: Store Antibody 1 separately at -80°C upon receipt. Avoid freeze-thaw cycles for 1% Digitonin.

Additional Materials

- Nuclease-free ddH₂O
- 80% ethanol
- BSA dry powder
- Vortexer
- Magnetic plate rack
- Centrifuge
- Orbital shaker
- Water bath or incubator
- Thermal cycler
- Analytical balance

Important Notes

- This kit is suitable for cultured cells or cells isolated from tissue samples. The recommended cell viability is > 85%.
- To avoid cell loss during wash processes, leave 3 ~ 5 uL liquid when removing the supernatant. A Horizontal centrifuge is easier for cell collection.

R-loop CUT&Tag Assay Protocol

The entire protocol takes 3 days to complete.

DAY 1

Prepare fresh buffers according to the tables below. Keep the buffers on ice.

10% BSA	
Component	Amount per sample
BSA dry power	50 mg
Nuclease-free ddH ₂ O	500 uL

Note: Store the freshly prepared BSA at 4°C. It will be used for buffer prep on Day 1 and Day 2.

Complete Wash Buffer	
Component	Amount per sample
Wash Buffer	100 uL
60 mM Spermidine	0.8 uL

Complete Antibody Buffer	
Component	Amount per sample
Antibody Buffer	83 uL
1% Digitonin	5 uL
1 % NP40	1 uL
60 mM Spermidine	0.8 uL
10% BSA	10 uL

1. Incubate Antibody 1

1.1. Count cells and check the cell viability (> 85%) by Trypan Blue staining. Aliquot 200,000 cells to a 1.5 mL tube per sample.

1.2. Centrifuge the cells at 600 x g at 4°C for 3 minutes. Remove and discard the supernatant.

1.3. Resuspend the cells in 100 uL Complete Wash Buffer. Centrifuge the cells at 600 x g at 4°C for 3 minutes. Remove and discard the supernatant.

1.4. Resuspend cells with 99 uL of Complete Antibody

Buffer and 1 uL of Antibody 1. Incubate the tube on an orbital shaker at 20 RPM overnight at 4°C.

DAY 2

Prepare fresh buffers according to the tables below. Keep the buffers on ice.

Complete Dig-150 Buffer	
Component	Amount per sample
Dig-150 Buffer (10x)	150 uL
1% Digitonin	75 uL
60 mM Spermidine	12 uL
10% BSA (from Day 1)	150 uL
1 % NP40	15 uL
Nuclease-free ddH ₂ O	1098 uL

Complete Dig-300 Buffer	
Component	Amount per sample
Dig-300 Buffer (10x)	110 uL
1% Digitonin	55 uL
60 mM Spermidine	9.2 uL
1 % NP40	11 uL
10% BSA (from Day 1)	110 uL
Nuclease-free ddH ₂ O	805 uL

2. Incubate Antibody 2

2.1. Centrifuge the cells at 600 x g for 5 minutes at 4°C. Remove and discard the supernatant. Avoid disturbing the cell pellet.

2.2. Resuspend the cells with 200 uL Complete Dig-150 Buffer by gentle pipetting. Incubate for 10 minutes at room temperature, gently mixing the cells once with a pipette at the 5-minute mark.

2.3. Repeat wash steps 2.1 and 2.2 twice (three washes in total).

2.4. Resuspend the cells with 99 uL Complete Dig-150 Buffer and 1 uL Antibody 2 by gentle pipetting. Transfer to a new 1.5 mL tube. Incubate the tube on an orbital shaker at 20 RPM for 1 hour at room temperature.

3. Incubate pA-Tn5 transposome

3.1. Repeat wash steps 2.1 - 2.2 three times in total.

3.2. Resuspend the cells with 98 uL Complete Dig-300 Buffer and 2 uL Assembled pA-Tn5 Transposome by gentle pipetting. Transfer to a new 1.5 mL tube. Incubate the tube on an orbital shaker at speed of 20 RPM for 1 hour at room

temperature.

1X	4°C	Hold
----	-----	------

4. Tagmentation of target DNA

4.1. Centrifuge cells at 600 x *g* for 5 minutes at 4°C. Remove and discard the supernatant. Avoid disturbing the cell pellet.

4.2. Resuspend the cells with 300 μ L Complete Dig-300 Buffer by gentle pipetting. Incubate for 10 minutes at room temperature, gently mixing the cells once with a pipette at the 5-minute mark.

4.3. Repeat wash steps 4.1 - 4.2 twice (three washes in total).

4.4. Resuspend the cells with 49.5 μ L Complete Dig-300 Buffer and 0.5 μ L Tagmentation Buffer by gentle pipetting. Transfer to a new 1.5 mL tube. Incubate at 37°C for 1 hour, with gentle mixing of the cells once using a pipette at the 30-minute time point.

5. Release of DNA

5.1. Mix Stop Buffer just before use. Add 6 μ L Stop Buffer and 1 μ L of Neutralize Buffer by gentle pipetting. Incubate at 55°C for 1 hour, gently mixing the cells with a pipette at the 20- and 40-minute time points.

5.2. Add 110 μ L DNA Beads and mix with a vortexer. Incubate at room temperature for 10 minutes.

5.3. Place the tube on a magnetic rack for 4 minutes to separate the beads. Remove and discard the supernatant.

5.4. Keep the tube on the magnetic rack, add 350 μ L of 80% ethanol without disturbing the beads. After 1 minute, remove and discard the supernatant. Repeat this step once for a total of two washes.

5.5. Air dry at 37°C for 3 min to remove the residual ethanol completely. Resuspend the beads with 20 μ L Nuclease-free ddH₂O. Transfer the beads to a new PCR tube.

6. Library amplification

6.1. Add 1 μ L P5 PCR Primer (20 μ M), 1 μ L P7 PCR Primer (20 μ M), and 22 μ L of PCR Mix (2x) to the PCR tube by gentle pipetting.

6.2. PCR Amplification using the following PCR program:

Cycle number	Temperature	Time
1X (initial heating)	72°C	5 minutes
1X (Denaturation)	98°C	30 seconds
15 ~ 17X	98°C	10 seconds
	60°C	20 seconds
	72°C	30 seconds
1X	72°C	5 minutes

DAY 3

7. Enrichment of library fragments

7.1. Resuspend and transfer the DNA beads to a new 1.5 mL tube. Add an extra 27 μ L DNA Beads to the PCR product. Mix thoroughly and incubate at room temperature for 10 minutes.

7.2. Place the sample on a magnetic rack for 4 minutes to separate the beads. Transfer the supernatant to a new 1.5 mL tube. Add 17 μ L DNA Beads to the supernatant. Mix and incubate at room temperature for 10 minutes.

7.3. Place the tube on a magnetic rack for 4 minutes to separate the beads. Remove and discard the supernatant.

7.4. Keep the tube on the magnetic rack, add 350 μ L of 80% ethanol without disturbing the beads. After 1 minute, remove and discard the supernatant. Repeat this step once for a total of two washes.

7.5. Allow the DNA beads to air dry completely. Add 10-20 μ L Nuclease-free ddH₂O to resuspend the beads and elute the DNA.

7.6. Use the magnetic rack to separate the beads. Transfer the supernatant (library) into a fresh tube.

8. Library quantification

8.1. Take 1 μ L library for DNA concentration measurement with a Qubit™ Fluorometer.

- If the total DNA amount is > 20 ng, the library is ready for sequencing.
- If the total DNA amount is < 20 ng, PCR amplify by 2-4 more cycles in step 6.

Kit performance

200,000 HEK293T cells were used to prepare R-loop CUT&Tag sequencing library using this kit according to the kit protocol. The R-loop profiling results were obtained as below (Fig.2, 3).




Fig.2. The distribution of R-loop peaks in the genome, showing characteristic enrichment at TSS.




Fig.3. The IGV visualization of R-loop CUT&Tag tracks on the AKT2 gene.

Troubleshooting

Problems	Possible reasons	Solutions
Low signal-to-background ratios in low-depth sequencing	Often caused by unhealthy cells or otherwise non-ideal input material.	1. Pre-treat cells with DNase or using flow cytometry to sort viable cells. 2. Sequence more deeply.
Low library amplification efficiency	Residual ethanol may inhibit the enzyme activity in step 5.3.	1. Increase the PCR cycle number. 2. Ensure ethanol is completely evaporated before DNA elution.

Appendixes

1. P5 and P7 PCR primer sequences (5'→3')

P5, Index 1	AATGATACGGCGACCACCGAGATCTACAC TATCCTCT TCGTCCGCGACGTC
P5, Index 2	AATGATACGGCGACCACCGAGATCTACAC GTAAGGA TCGTCCGCGACGTC
P5, Index 3	AATGATACGGCGACCACCGAGATCTACAC ACTGCATA TCGTCCGCGACGTC
P7, Index 1	CAAGCAGAAGACGGCATACGAGAT CTAGTACG GTCTCGTGGGCTCGG
P7, Index 2	CAAGCAGAAGACGGCATACGAGAT TTCTGCCT GTCTCGTGGGCTCGG
P7, Index 3	CAAGCAGAAGACGGCATACGAGAT GCTCAGGA GTCTCGTGGGCTCGG

Note: The sample multiplexing indexes are yellow highlighted and marked in green and red. A P5 PCR primer can be paired with a P7 PCR primer in any index combination for sample multiplexing in sequencer.

2. Library fragment schematic sequence:

5'-
AATGATACGGCGACCACCGAGATCTACAC **|||||**TCGTCCGCGACGTC
AGATGTGTATAAGAGACAG-**NNNNNN**-
CTGTCTCTTATACATCTCCGAGCCCACGAGAC **|||||**ATCTCGTATG
CCGTCTTCTGCTTG-3'

|||||: P5, Index, 8 bases
|||||: P7, Index, 8 bases
NNNNNN: Target DNA sequence