

Why Profile tRNA Modifications?

tRNAs are the most heavily modified RNA species in cells. These chemical modifications regulate tRNA structure, stability, decoding fidelity, and translational efficiency. Aberrant tRNA modification patterns are closely associated with cancer, neurological disorders, and metabolic diseases (Fig. 3).




Figure 3. tRNA modifications in diverse diseases.

Arraystar tRNA Modification-Seq

- Simultaneously identify the exact sites of m1A, m3C, m1G and m2,2G modifications

tRNAs modified with multiple methylations (m1A, m3C, m1G, m2,2G) can be simultaneously and quantitatively profiled at single-base resolution by tRNA Modification-Seq (m1A/m3C/m1G/m2,2G). These modifications disrupt Watson-Crick base pairing and can induce mismatch mutations due to misincorporation by reverse transcription detectable by sequencing.

To identify modification sites and quantify the methylation in a tRNA, the highly optimized tRNA-seq is performed on the tRNA with and without demethylase treatment. The modification is detected by methylation induced mutation index (MI) at each base position compared with its demethylase-treated tRNA. Among m1A, m3C, m1G, and m2,2G modifications, the m1G and m2,2G are further distinguished by their known distinct base positions in mammalian tRNAs (m1G at 9 and 37; m2,2G at 26 tRNA base positions).

Highlights

- Multiple methylation modification detection. Simultaneously profile m1A, m3C, m1G and m2,2G modifications in tRNAs
- High efficiency methylation removal. Efficiently remove internal and terminal modifications in the demethylase treated tRNA.
- Single-base resolution. Identify the exact methylation modification sites by comparing demethylase-untreated tRNA with demethylase-treated tRNA side by side.
- Comprehensive tRNA epitranscriptomic analysis. Analyze modifications sites, methylation levels, and tRNA expression all at once (Fig. 4).
- High performance tRNA sequencing. Highly optimized tRNA-seq with rigorous QC process, overcoming difficulties of regular sequencing for tRNA.
- Integrated tRNA Solutions. Work well together with our broad portfolio of tRNA research technologies: tRNA PCR array, LC-MS tRNA modification analysis, tRF&tiRNA-seq, and tRF&tiRNA PCR array.




Figure 4. Volcano plot showing differential tRNA methylation sites between two comparison groups.

References (will be a link to save space)



www.arraystar.com | info@arraystar.com | 888-416-6364